

PRODUCT REQUISITION FORM (PRF)

medi

Patient Name.....

I feel better
www.medi.de

Bed / Room No.

Below Knee (Calf)			mediven® thrombexin® 18			
Size	Extra Small	Small	Medium	Large	Large XWide	
Ankle	17 - 19	20 - 22	23 - 25	26 - 28	26 - 28	
Below Knee (Calf)	26 - 30	30 - 36	34 - 40	38 - 45	42 - 49	
Pl. Tick(✓) Size						




Circumference in cms.

Above Knee (Thigh)					mediven® thrombexin® 18				
Size	Extra Small	Small XWide	Medium	Medium XWide	Large	Large XWide	Extra Large	XLarge XWide	
Ankle	17-19	20-22	23-25	23-25	26-28	26-28	28-30	28-30	
Thigh	40-50	48-60	56-70	56-70	64-80	64-80	72-90	72-90	80-100
Pl. Tick(✓) Size									



Circumference in cms.

Below Knee (Calf)		mediven® struva® 23						
Size		I	II	III	IV	V	VI	VII
Ankle		17 - 19	19-21	21-23	23-25	25-27	27-29	29-31
Below Knee (Calf)		28-34	29-37	30-39	31-41	32-43	33-45	34-47
Pl. Tick(✓) Size								



Circumference in cms.

Above Knee (Thigh)		mediven® struva® 23						
Size		I	II	III	IV	V	VI	VII
Ankle		17-19	19-21	21-23	23-25	25-27	27-29	29-31
Below Knee (Calf)		28-34	29-37	30-39	31-41	32-43	33-45	34-47
Thigh		42-54	46-58	50-62	54-66	58-70	62-74	66-78
Pl. Tick(✓) Size								



Circumference in cms.

Below Knee (Calf)		mediven® struva® 35						
Size		I	II	III	IV	V	VI	VII
Ankle		17-19	19-21	21-23	23-25	25-27	27-29	29-31
Below Knee (Calf)		28-34	29-37	30-39	31-41	32-43	33-45	34-47
Pl. Tick(✓) Size								



Circumference in cms.

Above Knee (Thigh)		mediven® struva® 35						
Size		I	II	III	IV	V	VI	VII
Ankle		17-19	19-21	21-23	23-25	25-27	27-29	29-31
Below Knee (Calf)		28-34	29-37	30-39	31-41	32-43	33-45	34-47
Thigh		42-54	46-58	50-62	54-66	58-70	62-74	66-78
Pl. Tick(✓) Size								



Circumference in cms.


Measurement Method (All Leg Stocking)


- Please take and write the measurement as per the box on the right side having leg.
- Then check the size with the table with colour boxes on the left side and tick (✓) mark the required size.
- Proceed for order.





Phone : +91-33-2236-0368/67, Fax : +91-33- 2221-7335, e-mail : pushpanjali@vsnl.com
AN ISO 9001 : 2008 Certified Organisation • www.pushpanjaligroup.com

Available at :

Below Knee (Calf)			duomed®					Circumference in cms.
Size	Small	Medium	Large	Extra Large	Double Extra Large			
Ankle	19 - 21	22 - 24	25 - 27	28 - 30	30 - 32			
Below Knee (Calf)	28 - 34	32 - 38	36 - 42	40 - 46	42 - 50			
Pl. Tick(✓) Size								

Above Knee (Thigh)			duomed®					Circumference in cms.
Size	Small	Medium	Large	Extra Large	Double Extra Large			
Ankle	19 - 21	22 - 24	25 - 27	28 - 30	30 - 32			
Below Knee (Calf)	28 - 34	32 - 38	36 - 42	40 - 46	42 - 50			
Thigh	42 - 57	48 - 64	54 - 71	60 - 78	65 - 85			
Pl. Tick(✓) Size								

Below Knee (Calf)		mediven® Plus								Circumference in cms.
Size		I	II	III	IV	V	VI	VII		
Ankle		18-20	20-22	22-24	24-26	26-28	28-30	30-32		
Below Knee (Calf)		28-34	30-37	33-40	35-43	37-46	39-49	41-51		
Pl. Tick(✓) Size										

Above Knee (Thigh)		mediven® Plus								Circumference in cms.
Size		I	II	III	IV	V	VI	VII		
Ankle		18-20	20-22	22-24	24-26	26-28	28-30	30-32		
Below Knee (Calf)		28-34	30-37	33-40	35-43	37-46	39-49	41-51		
Thigh		43-57	45-62	49-67	53-72	56-77	60-81	64-85		
Pl. Tick(✓) Size										

mediven® armsleeves CCL 1* and 2- Lymphoedema Patients












TABLE - 2		armsleeve		armsleeve long		combined		combined long		handpiece
		710* <input type="checkbox"/> / 712 <input type="checkbox"/>		711* <input type="checkbox"/> / 713 <input type="checkbox"/>		730* <input type="checkbox"/> / 732 <input type="checkbox"/>		731* <input type="checkbox"/> / 733 <input type="checkbox"/>		720* <input type="checkbox"/> / 722 <input type="checkbox"/>
		armsleeve with shoulder cap and belt		armsleeve long with shoulder cap and belt		armsleeve with hand + shoulder cap + belt		armsleeve long with hand + shoulder cap + belt		glove/mittens
	738* <input type="checkbox"/> / 742 <input type="checkbox"/>		739* <input type="checkbox"/> / 743 <input type="checkbox"/>		740* <input type="checkbox"/> / 744 <input type="checkbox"/>		741* <input type="checkbox"/> / 745 <input type="checkbox"/>		723 <input type="checkbox"/>	

TABLE - 1		Circumference in cms.	Measurement for Sizes (in cms) for all the above styles						
			Size	II	III	IV	V	VI	VII
			cA	16-19	18-21	20-23	22-25	24-27	26-29
			cC	13-15	15-17	17-19	19-21	21-23	23-25
			cE	20-26	23-29	26-32	29-35	32-38	35-41
cG	22-29	25-32	29-36	33-40	37-44	41-48			

Box - 2 → Pl. Tick(✓) Size

Measurement Method (mediven® Armsleeve)

- First Take the measurement of cA, cC, cE & cG of the hand as per Table 1.
- Please tick (✓) the size in the box-2 as per table no. 1.
- Please tick (✓) in the box for the style as per Table no. 2
- Proceed for the order by matching the style as per Table no. 2 with the size marked (✓) in the box-2 of the table no. 1