

PRODUCT REQUISITION FORM (PRF)

Patient Name :

Bed / Room No. :


Doctor's Name :

medi

mediven® thrombexin® 18

Calf (Below Knee - AD)

Size	XS	S	M	L	LXW	XL	XXL* (on request)
Ankle cB	17 - 19	20 - 22	23 - 25	26 - 28	26 - 28	29 - 31	32 - 34
Calf cC	26 - 30	30 - 36	34 - 40	38 - 45	42 - 49	44 - 55	51 - 66
Tick (✓)	Left Right						




Circumference in cms.	
Left	Right

mediven® thrombexin® 18

Thigh (Above Knee - AG)

Size	XS	S	SXW	M	MXW	L	LXW	XL	XLW
Ankle cB	17-19	20-22	20-22	23-25	23-25	26-28	26-28	29-31	29-31
Thigh cG	40-50	48-60	56-70	56-70	64-80	64-80	72-90	72-90	80-100
Tick (✓)	Left Right								




Circumference in cms.	
Left	Right

mediven® struva® 23

Calf (Below Knee - AD)

Size	I	II	III	IV	V	VI	VII
Ankle cB	17 - 19	19-21	21-23	23-25	25-27	27-29	29-31
Calf cC	28-34	29-37	30-39	31-41	32-43	33-45	34-47
Tick (✓)	Left Right						




Circumference in cms.	
Left	Right

mediven® struva® 23

Thigh (Above Knee - AG)

Size	I	II	III	IV	V	VI	VII
Ankle cB	17-19	19-21	21-23	23-25	25-27	27-29	29-31
Calf cC	28-34	29-37	30-39	31-41	32-43	33-45	34-47
Thigh cG	42-54	46-58	50-62	54-66	58-70	62-74	66-78
Tick (✓)	Left Right						




Circumference in cms.	
Left	Right

mediven® struva® 35

Calf (Below Knee - AD)

Size	I	II	III	IV	V	VI	VII
Ankle cB	17-19	19-21	21-23	23-25	25-27	27-29	29-31
Calf cC	28-34	29-37	30-39	31-41	32-43	33-45	34-47
Tick (✓)	Left Right						




Circumference in cms.	
Left	Right

mediven® struva® 35

Thigh (Above Knee - AG)

Size	I	II	III	IV	V	VI	VII
Ankle cB	17-19	19-21	21-23	23-25	25-27	27-29	29-31
Calf cC	28-34	29-36	30-39	31-41	32-43	33-45	34-47
Thigh cG	42-54	46-58	50-62	54-66	58-70	62-74	66-78
Tick (✓)	Left Right						




Circumference in cms.	
Left	Right

duomed®

Calf (Below Knee - AD)

Size	S	M	L	XL	XXL
Ankle cB	19 - 21	22 - 24	25 - 27	28 - 30	30 - 32
Calf cC	28 - 34	32 - 38	36 - 42	40 - 46	42 - 50
Tick (✓)	Left Right				




Circumference in cms.	
Left	Right

duomed®



Thigh (Above Knee - AG)

Size	S	M	L	XL	XXL
Ankle cB	19 - 21	22 - 24	25 - 27	28 - 30	30 - 32
Calf cC	28 - 34	32 - 38	36 - 42	40 - 46	42 - 50
Thigh cG	42 - 57	48 - 64	54 - 71	60 - 78	65 - 85
Tick (✓)	Left Right				




Circumference in cms.	
Left	Right

...Contd/2

mediven® Plus				Calf (Below Knee - AD)					Circumference in cms.	
Size	I	II	III	IV	V	VI	VII		Left	Right
Ankle cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32			
Calf cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51			
Tick (✓)	Left									
	Right									
mediven® Plus				Thigh (Above Knee - AG)					Circumference in cms.	
Size	I	II	III	IV	V	VI	VII		Left	Right
Ankle cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32			
Calf cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51			
Thigh cG	43-48	45-62	49-66	53-60	56-64	60-68	64-72			
Tick (✓)	Left									
	Right									
Circaid® juxtalite® essentials (Lower Leg)				Calf (Below Knee - AD)					Circumference in cms.	
Size	S	M	M-X	L	L-X	XL	XL-X		XXL	Left
cC	26-36	33-44	44-54	41-51	51-61	37-48	48-61	53-64		
cB	19-22	22-27	25-32	41-51	25-32	30-37	30-37	35-42		
Tick (✓)	Left									
	Right									
Circaid® juxtalite® ankle foot wrap				Ankle (Below Knee - cA1)					Circumference in cms.	
Size	S			M					Left	Right
cA1	19-26			24-31						
Tick (✓)	Left									
	Right									
Circaid® juxtafit® essentials (Upper Leg)				Thigh (Above Knee - AG)					Circumference in cms.	
Size	XXS	M	M-X	L	L-X	XL			Left	Right
cG	40-52	46-62	54-70	60-80	60-90	77-101				
cE1	32-42	37-48	42-55	50-65	58-76	66-87				
cD	23-30	27-35	31-41	35-46	40-52	44-58				
Tick (✓)	Left									
	Right									
Circaid® juxtafit® essentials (Lower Leg)				Calf (Below Knee - AD)					Circumference in cms.	
Size	S	M	M-X	L	L-X	XL	XL-X		XXL	Left
cC	30-38	36-43	43-50	41-48	48-56	41-51	51-61	53-64		
cB1	20-28	25-33	33-40	30-38	38-46	33-41	41-51	43-53		
cB	15-23	20-28	20-28	25-33	25-33	30-38	30-38	33-41		
Tick (✓)	Left									
	Right									

mediven® harmony CCL 1* and 2-Lymphoedema Patients

TABLE - 1	TABLE - 2																																																								
	 •714H ☐	 716H ☐																																																							
	Circumference in cms. <table><tr><th>Left</th><th>Right</th></tr><tr><td> </td><td> </td></tr><tr><td> </td><td> </td></tr><tr><td> </td><td> </td></tr><tr><td> </td><td> </td></tr><tr><td> </td><td> </td></tr><tr><td> </td><td> </td></tr></table>	Left	Right													 •734H ☐																																									
		Left	Right																																																						
Box -2 → Tick(✓)	 •736H ☐																																																								
	 •738H ☐																																																								
	 742H ☐																																																								
	 •740H ☐																																																								
	 744H ☐																																																								
	 761H ☐																																																								
	Measurement for Sizes (in cms) for all the above styles																																																								
	<table><tr><th>Size</th><th>I</th><th>II</th><th>III</th><th>IV</th><th>V</th><th>VI</th><th>VII</th></tr><tr><td>cA</td><td>15.0-18.0</td><td>16.0-19.0</td><td>17.0-20.0</td><td>19.0-22.0</td><td>20.0-23.0</td><td>22.0-25.0</td><td>24.0-27.0</td></tr><tr><td>cC</td><td>13.5-14.5</td><td>14.5-15.5</td><td>16.0-17.5</td><td>17.5-19.0</td><td>19.0-20.5</td><td>20.5-22.5</td><td>22.5-24.5</td></tr><tr><td>cE</td><td>23.0-25.0</td><td>24.5-26.5</td><td>26.0-28.5</td><td>27.5-30.0</td><td>29.0-31.5</td><td>30.5-33.5</td><td>32.5-35.5</td></tr><tr><td>cG</td><td>25.5-28.0</td><td>28.0-30.5</td><td>30.0-33.0</td><td>32.0-35.0</td><td>33.5-36.5</td><td>36.0-39.5</td><td>38.5-41.5</td></tr><tr><td>Left</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>Right</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr></table>	Size	I	II	III	IV	V	VI	VII	cA	15.0-18.0	16.0-19.0	17.0-20.0	19.0-22.0	20.0-23.0	22.0-25.0	24.0-27.0	cC	13.5-14.5	14.5-15.5	16.0-17.5	17.5-19.0	19.0-20.5	20.5-22.5	22.5-24.5	cE	23.0-25.0	24.5-26.5	26.0-28.5	27.5-30.0	29.0-31.5	30.5-33.5	32.5-35.5	cG	25.5-28.0	28.0-30.5	30.0-33.0	32.0-35.0	33.5-36.5	36.0-39.5	38.5-41.5	Left								Right							
Size	I	II	III	IV	V	VI	VII																																																		
cA	15.0-18.0	16.0-19.0	17.0-20.0	19.0-22.0	20.0-23.0	22.0-25.0	24.0-27.0																																																		
cC	13.5-14.5	14.5-15.5	16.0-17.5	17.5-19.0	19.0-20.5	20.5-22.5	22.5-24.5																																																		
cE	23.0-25.0	24.5-26.5	26.0-28.5	27.5-30.0	29.0-31.5	30.5-33.5	32.5-35.5																																																		
cG	25.5-28.0	28.0-30.5	30.0-33.0	32.0-35.0	33.5-36.5	36.0-39.5	38.5-41.5																																																		
Left																																																									
Right																																																									

Box -2 → Tick (✓)

Measurement Method (For all Stockings range)

- Please mention the measurement as per the box on the right side showing the leg
- Check the size with the table matching the colour boxes on the left side and tick (✓) the required size
- Proceed with the order

Measurement Method (mediven harmony®)

- First mention the measurement of cA, cC, cE & cG of the hand as indicated in TABLE 1
- Please tick (✓) the size in the box-2 as per TABLE 1
- Please tick (✓) in the box for the style as per TABLE 2
- Proceed with the order by matching the style as per TABLE 2 with the size marked (✓) in the box 2 of TABLE 1

Patient Care Centre : **PUSHPANJALI** Hi-tech Rehab Centre, 17, Ganesh Chandra Avenue, Kolkata-700 013



PUSHPANJALI medi India Pvt Ltd



MEDICAL COMPRESSION | ORTHOPEDICS | PHYSIOTHERAPY | REHABILITATION
16, Ganesh Chandra Avenue, Gandhi House, Kolkata-700 013, West Bengal, INDIA
☎ : +91-33-4040 1300 📠 : 91633 60368 🌐 : www.pushpanjaligroup.com ✉ : info@pushpanjaligroup.com 📱 : pushpanjali medi india

Available at :