## Spinomed®V/Spinomed®II



Back orthosis for vertebral extension in osteoporosis and strengthening of the core musculature

## **QUALITY FEATURES**



Effortless donning and doffing thanks to pre-shaped shoulder

Back splint can be cold-formed without tools for individual adaptation to the spine



straightening and pain relief in accordance with therapeutic needs.



Swivel joints on the abdominal fastener allows optimised fit.



Pivot joint for optimum freedom of movement, secure fit and high



Tool-free fastening of the strap eyelets.




Bio-feedback mechanism for postural correction



Spinomed ® V

## **QUICK APPLICATION**





- Strengthen the back and abdominal muscles within six months (by up to 72% and 56%)
- relieve pain (by up to 47%).
- reduce the kyphosis angle (by up to 11%).
- reduce the tendency to sway (by up to 23%).

Ref. : Friedrich-Alexander University, Erlangen-Nurer randomised, controlled and partially blind study (Kemmler, W : Study on th Spinomed Back Brace with regard to Swaying and Compliance compared witl Competitors', 2014, unpublised). One of the study's most important endpoint

## Spinomed® II



Scientifically proven effectiveness - The evidence based study of Spinomed® was published in the "American Journal of Physical Medicine & Rehabilitation", edition March 2004., 83:177-186 Authors: Michael Pfeifer, MD; Bettina Begerow, PhD; Helmut W Minne, MD



	Measure	XS	S	M	L	XL*	XXL*
back length (cm)	bL	up to 42	41 - 48	47- 54	53 - 60	59 - 66	65 - 72

Measurement: Length from below C7 to sacrum along spinal curvature





PUSHPANJALI medi India Pvt Ltd MEDICAL COMPRESSION | ORTHOPEDICS | PHYSIOTHERAPY | REHABILITATION | We see a collection of the control of the



