

# Hot & Cold Therapy

## Polar Frost® Cold Gel from Finland

Made specifically for Mettler Electronics, Polar Frost Cold Gel is specially formulated to provide soothing relief for soft tissue injuries, stress injuries, muscle tension, inflammation and stiffness. Aloe Vera has an anti-inflammatory effect and it keeps the skin smooth and moist throughout treatment. Polar Frost is available in a variety of sizes that includes a roll-on bottle for fingers-free application. Massage a thin layer of Polar Frost Cold Gel into the painful area. Repeat application every three hours for 3-5 days, as needed. This gel absorbs fast and works quickly. Try Polar Frost and see what the other products are lacking.



Polar Frost is manufactured under a strictly standardized production process using only the highest quality materials.

### Indications

Provides advanced relief for :

- Arthritis pain • Joint pain • Back pain • Knee pain • Bunion pain
- Muscle pain • Foot pain • Neck pain • Heel pain • Shoulder pain

Item#	Description	Size
1881	Polar Frost, pump bottle	500 ml. (17 oz.)
1882	Polar Frost, tube	150 ml. (5 oz.)
1883	Polar Frost, roll-on	75 ml. (2.5 oz.)
1884	Polar Frost	1 gallon
1885	Polar Frost	4 x 1 gallon
1886	Polar Frost sample packages	200 pieces

## Electrodes

### EZ Trode® Self-Adhesive Electrodes (Now with Aloe Vera)

EZ Trode self-adhesive electrodes with Aloe Vera have a thin silver layer that provides superior conductivity. They have soft-rubber connectors, are long-lasting and come in a special resealable pouch that preserves freshness while being easy to open and close. EZ Trodes come in four convenient sizes.

#### Item# Description

- 2221 EZ Trode – 2" diameter round electrodes with lead wires, case of ten packages (four electrodes/pkg.)
- 2222 EZ Trode – 2.75" diameter round electrodes with lead wires, case of ten packages (four electrodes/pkg.)
- 2223 EZ Trode – 2" x 5" rectangular electrodes with lead wires, case of 10 packages (two electrodes/pkg.)
- 2224 EZ Trode – 2" square electrodes with lead wires, case of 10 packages (four electrodes/pkg.)



### V Trode™ Self-Adhesive Electrodes

V Trode self-adhesive neurostimulation electrodes feature soft rubber lead wire attachments that make it easy to attach them to the electrode lead wires. A long lasting gel with a thin carbon film layer provides strong adhesion to the skin. V Trodes come in four convenient sizes.

#### Item# Description

- 2702 V Trode – 2" diameter round electrodes with lead wires, case of 10 packages (four electrodes/pkg.)
- 2703 V Trode – 2.75" diameter round electrodes with lead wires, case of 10 packages (four electrodes/pkg.)
- 2704 V Trode – 2" x 4" oval electrodes with lead wires, case of 10 packages (four electrodes/pkg.)
- 2705 V Trode – 2" square electrodes with lead wires, case of 10 packages (four electrodes/pkg.)



# Ultrasound Accessories

## Ultrasound Applicators

### Item# Description

- 7310 Sonicator 716, 730, 992 or 994 applicator (10 cm<sup>2</sup>/ 1 MHz)
- 7335 Sonicator 730 applicator (5 cm<sup>2</sup>/ 3.3 MHz)
- 7305 Sonicator 715 or 730 applicator (5 cm<sup>2</sup>/ 1 MHz)
- 7513 Sonicator Plus 930, 992 or 994 applicator (5 cm<sup>2</sup> / 1.0 & 3.2 MHz)
- 7331 Sonicator 730, 992 or 994 pencil style applicator (1 cm<sup>2</sup>/ 3.3 MHz)

New applicators for the Sonicator 740 with the bright blue coupling indicator light:

- 7413 5 cm<sup>2</sup> / 1.0 & 3.2 MHz
- 7410 10 cm<sup>2</sup> / 1 MHz
- 7431 1 cm<sup>2</sup> / 3.3 MHz



## Sonicator Plus 920 and 940 Applicators

### Item# Description

- 9401 5.5 cm<sup>2</sup>, 1/3 MHz applicator for the Sonicator Plus 940
- 9402 0.9 cm<sup>2</sup>, 1/3 MHz applicator for the Sonicator Plus 940
- 9201 5.5 cm<sup>2</sup>, 1/3 MHz applicator for the Sonicator Plus 920
- 9202 0.9 cm<sup>2</sup>, 1/3 MHz applicator for the Sonicator Plus 920

## Stimulator Accessories

### Pencil Electrode Set

The Pencil Electrode is used for manual stimulation procedures where pin-point control is required. It comes in its own carry case with four stainless steel tips from the size of a quarter down to a point. The Pencil Electrode has an On /Off switch so that stimulation can be turned on and off. It requires a banana type connector to be plugged into the electrode cable set. It is compatible with all Mettler stimulators.

### Item# Description

- 2023 Pencil electrode set with push button stimulation control, (includes handle, 4 different sizes of stainless steel spot electrode tips, and carrying case)



### Sponge Inserts for Sponge Electrodes

- 2001 24 Sponge inserts (2" x 2")
- 2003 24 Sponge inserts (4" x 4")
- 2005 12 Sponge inserts (3.5" x 7")
- 2007 12 Sponge inserts (8" x 10")

### Sponge Electrodes

- 2000 4 Sponge electrodes (2" x 2")
- 2002 4 Sponge electrodes (4" x 4")
- 2004 1 Sponge electrode (3.5" x 7")
- 2006 1 Sponge electrode (8" x 10")

### Electrode Straps

- 2008 4 Electrode straps (24")
- 2009 4 Electrode straps (48")

