

CPM

CPM 8000

Mettler Continuous Passive Motion

Introducing the Mettler CPM 8000. It provides complete Automatic motion of the knee. It is ideal for all kinds of knee CPM patients including orthopedic, pediatric and bariatric.

Features

- Easy to use • Adjustable • Patient lock-out
- Wide carriage for all sizes of patients
- 2 Year Warranty • Programmable
- Reversible footplate • Easy to use LCD control unit



SPECIFICATIONS

Input :	100 to 230 VAC , 50/60 Hz
Weight :	33 lbs., (15 kg)
Dimensions :	47" x 12" x 16.5", (120 x 30 x 42 (cm))
Treatment Time :	Continuous
Speeds	
Min. (100%)	70 s ± 10% (- 10° to 120°) 75 s ± 10% (120° to - 10°)
Mean (60%)	120 s ± 10% (- 10° to 120°) 120 s ± 10% (120° to - 10°)
Max. (30%)	240 s ± 10% (- 10° to 120°) 250 s ± 10% (120° to - 10°)
Working Angles :	- 10° (± 2°) to 120° (± 2°)
Range of Motion :	For H = 8 from -10° to 120°
(ROM)	For H = 5 from -10° to 115°
	For H = 0 from -10° to 110°

Hot & Cold Therapy

ThermalSoft® Gel Hot and Cold Packs

ThermalSoft Gel packs were designed with one thing in mind. Quality! Mettler set out to find the strongest most durable hot/cold pack on the market. Our innovation with our Double Wall construction did just that! This allows the packs to be softer and more durable than most other products. Unlike other cold packs on the market, ThermalSoft Gel Packs can do both hot and cold from the same pack. Heat them in the microwave for heat therapy or for cold therapy put them in the freezer. Due to the technology inside the ThermalSoft Gel Packs stay soft, even at freezer temperatures and retain their hot or cold therapeutic benefits for extended periods of time. The large gel pieces mold to your body to deliver soothing heat or cold therapy. All ThermalSoft Gel packs carry a one-year warranty.



Item#	Description	Size
9831	Standard Gel Pack	11" x 14"
9832	Half Size Gel Pack	7.5" x 11"
9833	Cervical Gel Pack	23" x 8" contour
9834	Extra Large Gel Pack	11" x 21"

ThermalSoft® Durapak Hot and Cold Packs

ThermalSoft Durapak hot and cold packs maintain an even temperature over the entire pad. Their soft gel contents are 100% leak-proof and contain no liquids or dyes. ThermalSoft Durapak hot and cold packs are biodegradable and non-caustic. And, because ThermalSoft Durapaks are reusable, they save you money, too. For heat therapy, microwave or heat the pack in water. ThermalSoft Durapaks keep their pad like shape even at high temperatures and will not liquefy when heated. When used for cold therapy, ThermalSoft Durapaks stay flexible when frozen. They bend to fit around treatment areas and maintain their thickness no matter how severe the bend.



Item#	Description	Size
9841	Small	4" x 6"
9842	Standard	5" x 10"
9843	Large	8" x 11"
9844	Extra Large	12" x 15"



Hot & Cold Therapy

Polar Frost® Cold Gel from Finland

Made specifically for Mettler Electronics, Polar Frost Cold Gel is specially formulated to provide soothing relief for soft tissue injuries, stress injuries, muscle tension, inflammation and stiffness. Aloe Vera has an anti-inflammatory effect and it keeps the skin smooth and moist throughout treatment. Polar Frost is available in a variety of sizes that includes a roll-on bottle for fingers-free application. Massage a thin layer of Polar Frost Cold Gel into the painful area. Repeat application every three hours for 3-5 days, as needed. This gel absorbs fast and works quickly. Try Polar Frost and see what the other products are lacking.



Polar Frost is manufactured under a strictly standardized production process using only the highest quality materials.

Item#	Description	Size
1881	Polar Frost, pump bottle	500 ml. (17 oz.)
1882	Polar Frost, tube	150 ml. (5 oz.)
1883	Polar Frost, roll-on	75 ml. (2.5 oz.)
1884	Polar Frost	1 gallon
1885	Polar Frost	4 x 1 gallon
1886	Polar Frost sample packages	200 pieces

Indications

Provides advanced relief for :

- Arthritis pain
- Joint pain
- Back pain
- Knee pain
- Bunion pain
- Muscle pain
- Foot pain
- Neck pain
- Heel pain
- Shoulder pain

Electrodes

EZ Trode® Self-Adhesive Electrodes (Now with Aloe Vera)

EZ Trode self-adhesive electrodes with Aloe Vera have a thin silver layer that provides superior conductivity. They have soft-rubber connectors, are long-lasting and come in a special resealable pouch that preserves freshness while being easy to open and close. EZ Trodes come in four convenient sizes.

Item# Description

- 2221 EZ Trode – 2" diameter round electrodes with lead wires, case of ten packages (four electrodes/pkg.)
- 2222 EZ Trode – 2.75" diameter round electrodes with lead wires, case of ten packages (four electrodes/pkg.)
- 2223 EZ Trode – 2" x 5" rectangular electrodes with lead wires, case of 10 packages (two electrodes/pkg.)
- 2224 EZ Trode – 2" square electrodes with lead wires, case of 10 packages (four electrodes/pkg.)



V Trode™ Self-Adhesive Electrodes

V Trode self-adhesive neurostimulation electrodes feature soft rubber lead wire attachments that make it easy to attach them to the electrode lead wires. A long lasting gel with a thin carbon film layer provides strong adhesion to the skin. V Trodes come in four convenient sizes.

Item# Description

- 2702 V Trode – 2" diameter round electrodes with lead wires, case of 10 packages (four electrodes/pkg.)
- 2703 V Trode – 2.75" diameter round electrodes with lead wires, case of 10 packages (four electrodes/pkg.)
- 2704 V Trode – 2" x 4" oval electrodes with lead wires, case of 10 packages (four electrodes/pkg.)
- 2705 V Trode – 2" square electrodes with lead wires, case of 10 packages (four electrodes/pkg.)

