Medical Compression Garments medi

For realizing successful surgical outcomes after plastic-cosmetic surgery

I feel better

lipomed[®] basic





Quality Features



Different Models

Measurement

Circumference		XS	S	M	L	XĽ
below breast	(cN)	59 - 72	69 - 82	79 - 92	89 - 102	99 - 112
waist	(cM)	57 - 69	67 - 79	77 - 89	87 - 99	97 - 109
hip	(cl)	84 - 100	92 - 108	100 - 116	108 - 124	116 - 132
thigh	(cG)	47 - 57	54 - 64	61 - 71	68 - 78	75 - 85
mid-thigh	(cF)	38 - 46	42 - 50	46 - 54	50 - 58	54 - 62
below knee	(cD)	27 - 34	30 - 37	34 - 41	38 - 45	42 - 49
ankle	(cB)	17 - 19	19 - 22	22 - 26	26 - 30	

lipomed bra





Baumwolle-Special material that self adapts the body contour



Hook-and-eye closures in the straps



Combination front closure

Circumference in cm	S	M	L	XL
Below breast	70-80	80-90	90-105	105 - 125

Quality

Features



Colour coding system Markings on the Velcro fasteners make it easier to put on



lipomed[®] face

Measuring chart lipomed face					
Circumference in cm	S	L			
Neck circumference	< 40	>38			