

The economical solution for venous problems

QUALITY FEATURES




CCL 1/2


Compression Class



Clima Fresh



Perfect Fit



Graduated Compression

QUALITY CERTIFICATES




CE mark



Quality label (RAL)



Textile trust
Quality seal




HOSY-
Hohenstein
Testing Device



DIN EN ISO 9001:2008
DIN EN 13485


AVAILABLE MODELS




Calf

Thigh

Panty




Close toe



Open toe


MEASUREMENT

Sizes	cB Ankle measurement	cC Calf measurement	cG Thigh measurement	measurements for silicon topband
S	19 - 21 cm	28 - 34 cm	42 - 57 cm	47 - 54 cm
M	22 - 24 cm	32 - 38 cm	48 - 64 cm	52 - 60 cm
L	25 - 27 cm	36 - 42 cm	54 - 71 cm	56 - 65 cm
XL	28 - 30 cm	40 - 46 cm	60 - 78 cm	62 - 72 cm
XXL	30 - 32 cm	42 - 50 cm	65 - 85 cm	65 - 76 cm



Correct graduation of Pressure depends on correct selection of the ankle size

(I) **Emter. M.** Compression Stocking management of patients suffering from Venous Disorder (II) Phlebology Net(5 Abstracts) (Stockings)




Varicose Veins

weakened walls

valve does not
close properly

inflammation



Healthy Veins

blood flow

valve closes



PUSHPANJALI medi India Pvt Ltd

MEDICAL COMPRESSION | ORTHOPEDICS | PHYSIOTHERAPY | REHABILITATION



We serve all!
An ISO 9001:2015 Company



SCAN FOR MORE DETAILS