


mediven® thrombexin® 18

medi

Medical thrombosis prophylaxis stocking with 18 mmHg at the ankle



QUALITY FEATURES



Graduated Compression




Silicone beaded
top band



Anti-bacterial Toe



3D Fleecy heel



Perfect size indicator



Compression strength

INTELLIGENT TOP BAND

QUALITY CERTIFICATES

CE mark



Quality label
(RAL)



Textile trust
Quality seal



DIN EN ISO:9001:2008
DIN EN 13485



HOSY-
Hohenstein
Testing Device



3 DIFFERENT MODELS



MEASUREMENT

Knee (AD)

Size	Leg measurements* (circumference in cm)			Colour code
Calf (AD)	Ankle cB	Calf cC	Heel	Thigh
XS	17 – 19	26 – 32	orange	–
S	20 – 22	30 – 38	green	–
M	23 – 25	34 – 44	violet	–
L	26 – 28	38 – 48	blue	–
L-X	26 – 28	44 – 55	grey	–
XL	29 – 31	44 – 55	red	–
XXL*	32 – 34	51 – 66	brown	–


Thigh (AG)

Size	Ankle cB	Thigh cG	Colour code	Colour code
XS	17 – 19	40 – 50	orange	orange
S	20 – 22	48 – 60	green	green
S-X	20 – 22	56 – 70	green	grey
M	23 – 25	56 – 70	violet	violet
M-X	23 – 25	64 – 80	violet	grey
L	26 – 28	64 – 80	blue	blue
L-X	26 – 28	72 – 90	blue	grey
XL	29 – 31	72 – 90	red	red
XL-X	29 – 31	80 – 100	red	grey

Hip (AG)

Size	Ankle cB	Thigh cG	Colour code	Colour code
S	20 – 22	48 – 60	green	green
M	23 – 25	56 – 70	violet	violet
L	26 – 28	64 – 80	blue	blue
L-X	26 – 28	72 – 90	blue	grey

*on request



(I) **Stephen Thomas et al**, Graduated external compression and prevention of Deep Vein Thrombosis

(II) **Jennifer Hill, 2007**, Reducing the risk of Venous Thrombolism (Deep Vein Thrombosis & Pulmonary Embolism)



PUSHPANJALI medi India Pvt Ltd

MEDICAL COMPRESSION | ORTHOPEDICS | PHYSIOTHERAPY | REHABILITATION



SCAN FOR MORE DETAILS

+91-33-4040 1300 | 91633 60368 | www.pushpanjaligroup.com | info@pushpanjaligroup.com | pushpanjali medi india | shop-pushpanjali.com

*All prices are subject to change without notice, and are not guaranteed. Please check during purchase / sale.