# circaid<sup>®</sup> juxtafit<sup>®</sup> essentials arm Adjustable compression device for the arm





# **QUALITY FEATURES**







Adjustable, measurable and re-adjustable compression

**Built-In Pressure System (BPS)** 

# ADJUSTABLE COMPRESSION

10 - 20 mmHg | 20 - 30 mmHg | 30 - 40 mmHg

## **MATERIAL COMPONENTS:**

53% Nylon, 40% Polyurethane, 7% Elastane

# **QUALITY CERTIFICATES**









#### **BENEFITS & FEATURES:**

- Compression therapy essential part of the complex lymphoedema therapy
- Easy handling simple donning and doffing allows for selfmanagement of the patient
- Compression solution for optimising decongestion results and stabilising oedema fluctuations



## **CIRCAID IUXTAFIT ESSENTIALS ARM KIT INCLUDES:**

- 1 adjustable compression garment
- 1 circaid Built-in-Pressure system (BPS) guide
- 1 elbow band
- 1 circaid cover up arm

\	Circumferences in cm	1	II	III	IV	V	VI
7	cG	28 – 32	30 - 34	32 - 37	35 – 40	38 - 44	42 – 48
	сE	25 – 28	27 – 30	29 – 32	31 – 34	33 – 36	36 – 40
	сС	15 – 17	17 – 19	19 – 21	21 – 23	23 – 25	25 – 28

### medi tip:

Due to product features we recommend circaid kit for :

Men / women

- Patients who have challenges with medical compression garments
- Up to 24 hour compression therapy
- Mobile / immobile patients

	Length in cm	short	long
	ℓC-E	22 – 25	> 25
	Actual garment length in cm	short	long
	ℓC – G	43	48
	Chudian and I have the		1 1 1 1



cB	Circumferences in cm	xs	S	M	L	XL
	cB Palm	14 – 17	17 – 20	20 – 23	23 – 26	26 – 29

Studies: 1. Lymphoedema Framework. Best Practice for the Management of Lymphoedema. International consensus. London: MEP Ltd, 2006.

2. Protz K et al. Comparison of handling, pressure and wearing comfort of different compression edvices for decongestion therapy. Veins and Lymphatics 2018.



Y/AL/ medi India Pvt Ltd MEDICAL COMPRESSION | ORTHOPEDICS | PHYSIOTHERAPY | REHABILITATION





🕓 : +91-33-4040 1300 🖥 🖫 : 91633 60368 🧼 : www.pushpanjaligroup.com 🖂 : info@pushpanjaligroup.com 🖫 : pushpanjali medi india 📜 : shop-pushpanjali.com